
BEC Monitor Documentation

Release 0.0.1

Zachary Glassman

June 04, 2015

CONTENTS

1	Usage	3
1.1	Fitobject	3
1.2	Datatablewidget	7
1.3	Auxfuncwidget	7
1.4	Ipython	8
1.5	Auxwidgets	9
1.6	Subroutines	9
1.7	Dataplots	9
1.8	Fitmodels	10
1.9	Image	10
1.10	SpinorMonitor	11
1.11	Optionswidgets	13
1.12	Auxfunctions	14
1.13	Visualplotterwidget	14
2	Indices and tables	17
	Index	19

This software is used in the Lett Lab of the Laser Cooling and Trapping group at NIST/JQI for experiments on Spinor Na Bose-Einstein condensates.

Install the required packages

- numpy
- scipy
- pyqtgraph
- pyqt
- lmfit
- pandas

To use just run SpinorMonitor.py

Contents:

1.1 Fitobject

Contents:

1.1.1 fit_object

class Fitobject.**fit_object** (*index, params, type_of_fit, roi, data*)
fit object holds all the information for a single fit_sequence

Parameters

- **index** (*string*) – Shot number
- **params** (*dictionary*) – dictionary of Parameters objects containing fit parameters
- **type_of_fit** (*string*) – type of fit to be performed
- **roi** (*list*) – region of interest to crop data for fit
- **data** (*numpy array*) – numpy array of image to be analyzed

Methods

BEC_num (*scalex*, *scaley*)

get number of BEC atoms from fit from equation

$$N = \left(\frac{2\pi}{3\lambda^2} \right) \frac{2\pi A}{5} R_x R_y$$

Parameters

- **scalex** – x scale of pixel
- **scaley** – y scale of pixel

Variables

- **A** – fitted Thomas-Fermi amplitude
- **Rx** – fitted Thomas-Fermi x radius
- **Ry** – fitted Thomas-Fermi y radius
- **sigma** – optical density

Returns atom number

BEC_num_1 (*scalex*, *scaley*, *A*, *dx*, *dy*)

helper function for BEC num get number of BEC atoms from fit from equation

$$N = \left(\frac{2\pi}{3\lambda^2} \right) \frac{2\pi A}{5} R_x R_y$$

Parameters

- **scalex** – x scale of pixel
- **scaley** – y scale of pixel
- **A** – fitted Thomas-Fermi amplitude
- **dx** – fitted Thomas-Fermi x radius
- **dy** – fitted Thomas-Fermi y radius

Variables **sigma** – optical density

Returns atom number

TF_2D ()

two dimensional Thomas Fermi which is not normalized of the form:

$$TF = A \max \left\{ \left[1 - \left(\frac{x_c}{dx} \right)^2 - \left(\frac{y_c}{dy} \right)^2 \right], 0 \right\}^{3/2}$$

Variables

- **x0** – absolute x center
- **y0** – absolute y center
- **xc** – rotated x center
- **yc** – rotated y center
- **theta** – angle relative to x axis

- **A** – amplitude
- **dx** – Thomas-Fermi radius on rotated x axis
- **dy** – Thomas-Fermi radius on rotated y axis
- **off** – offset

Therm_num (*scalex, scaley*)

get number of BEC atoms from fit from equation

$$N = \left(\frac{2\pi}{3\lambda^2} \right) \frac{2\pi A}{5} R_x R_y$$

Parameters

- **scalex** – x scale of pixel
- **scaley** – y scale of pixel

Variables

- **A** – fitted Gaussian amplitude
- **Rx** – fitted Gaussian x standard deviation
- **Ry** – fitted Gaussian y standard deviation
- **sigma** – optical density

Returns atom number

bimod2min (*params*)

function to minimize, need to subtract offset since included in both terms

create_vecs (*roi*)

create vectors scaled by pixel size

Parameters **roi** (*list*) – region of interest list

Return X x vector from meshgrid

Return Y y vector from meshgrid

fit_image ()

fit corrected image with parameters from params

gauss_2D ()

two dimensional Gaussian which is not normalized of the form:

$$G = A \exp \left(-\frac{(x - x_c)^2}{2dx^2} - \frac{(y - y_c)^2}{2dy^2} \right) + Off$$

Variables

- **x0** – absolute x center
- **y0** – absolute y center
- **xc** – rotated x center
- **yc** – rotated y center
- **theta** – angle relative to x axis
- **A** – amplitude
- **dx** – standard deviation on rotated x axis

- **dy** – standard deviation on rotated y axis
- **off** – offset

get_angled_line (*x0, y0, theta*)

get angled line for angle theta with formulas

$$x_c = (x - x_0) \cos(\theta) - (y - y_0) \sin(\theta)$$

$$y_c = (x - x_0) \sin(\theta) - (y - y_0) \cos(\theta)$$

Parameters

- **x0** – absolute x center
- **y0** – absolute y center
- **xc** – rotated x center
- **yc** – rotated y center
- **theta** – angle relative to x axis

line_profile ()

calculate line profile, with zeroes to make full image

Returns two-dimensional array which has padding outside of the region of

interest and can be summed for profiles.

multiple_fits ()

function to fit sequentially with input defined from SpinorMonitor we may need to take parameters of previous fit!! do fit, update values, do next fit

partial_TF_2D (*xc, yc, A, dx, dy*)

two dimensional non-rotated Thomas Fermi which is not normalized of the form:

$$TF = A \max \left\{ \left[1 - \left(\frac{x_c}{dx} \right)^2 - \left(\frac{y_c}{dy} \right)^2 \right], 0 \right\}^{3/2}$$

Parameters

- **xc** – absolute x center
- **yc** – absolute y center
- **A** – amplitude
- **dx** – Thomas-Fermi radius on rotated x axis
- **dy** – Thomas-Fermi radius on rotated y axis
- **off** – offset

process_results (*scalex, scaley*)

process results of fit and allow output return dictionary scale with the appropriate pixel values after fit

sg2min (*params*)

stern gerlach function to minimize

stern_gerlach_2D ()

2 dimensional three thomas fermi distributions

subtract_background ()

Subtract background from image looking at first and last 20 rows of the initial image far away from experiment

1.2 Datatablewidget

Contents:

1.2.1 DataTable

class Datatablewidget.**DataTable** (*parent=None*)
tabbed tables to show system parameters and fitted parameters

Methods

update_pandas_table (*df*)
update tables, check if cols are different

1.3 Auxfuncwidget

Contents:

1.3.1 AuxillaryFunctionContainerWidget

class Auxfuncwidget.**AuxillaryFunctionContainerWidget** (*parent=None*)
class for displaying container of auxillary function widgets will hold a stacked layout of all auxillary functions

Methods

add_element (*name*)
convenicne function to create function widget and add to proper dictionaries
re_import ()

1.3.2 AuxillaryFunctionWidget

class Auxfuncwidget.**AuxillaryFunctionWidget** (*func, parent=None*)
class holding function and entry information

Methods

calculate ()
generate_info_widgets ()
generate info sublayouts
generate_params_widgets ()
generate parameter sublayout and return layout
get_params ()

1.4 Ipython

Contents:

1.4.1 PlotObj

class `Ipython.PlotObj`
class to hold SpinorPlot objects

Methods

add_plot (*plot, name*)
Add plot to dictionary of plots to update

update (*var_dict*)
update all plots in dictionary

1.4.2 QIPythonWidget

class `Ipython.QIPythonWidget` (*customBanner=None, *args, **kwargs*)
Convenience class for a live IPython console widget. This widget lives within the main GUI

Attributes

<code>custom_control</code>	
<code>custom_page_control</code>	

Methods

clearTerminal ()
Clears the terminal

executeCommand (*command*)
Execute a command in the frame of the console widget

printText (*text*)
Prints some plain text to the console

pushVariables (*variableDict*)
Given a dictionary containing name / value pairs, push those variables to the IPython console widget

1.4.3 QIPythonWidgetContainer

class `Ipython.QIPythonWidgetContainer` (*parent=None*)
Ipython container class for multi-threading

Methods

1.4.4 SpinorPlot

class `Ipython.SpinorPlot` (*func, name=None, xaxis=None, yaxis=None*)
class to plot with updating stuff

Methods

get_vars (*var_dict*)

set_axis ()

update_plot (*var_dict*)

1.5 Auxwidgets

Contents:

1.5.1 FingerTabBarWidget

class `Auxwidgets.FingerTabBarWidget` (*parent=None, *args, **kwargs*)
Class to implement tabbed browsing for options

Methods

paintEvent (*event*)

tabSizeHint (*index*)

1.5.2 TextBox

class `Auxwidgets.TextBox`
custom textbox, mostly QTextEdit, with some added functions

Methods

output (*x*)

1.6 Subroutines

Contents:

1.7 Dataplots

Contents:

1.7.1 ImageWindow

class `Dataplots.ImageWindow` (*parent=None*)
Image View with custom ROI

Attributes

<code>lastFileDir</code>	
--------------------------	--

Methods

add_lines (*results*)
add lines to plot, input it numpy array which is then summed

popup ()
function to start popup window object

setImage (*im*)
set image

updatePlot ()
updates plot, can only be called once plot is initialized with image

1.7.2 PlotGrid

class `Dataplots.PlotGrid` (*parent=None*)

Methods

1.8 Fitmodels

Contents:

1.9 Image

Contents:

1.9.1 IncomingImage

class `Image.IncomingImage`
check for images, if found obtain image and send back to main GUI

Methods

newImage ()
This function checks in directory for new image with proper name if found, it reads it in and then deletes it

run()

every second search folder for new images, if found get image and emit back to main gui for processing

1.9.2 ProcessImage

class Image.**ProcessImage** (*data, exp_data, options, path, run*)

Processing object for threading purposes @parameters

data: numpy array options: list of options for fit parameters

[params,

type_of_fit, ROI, index

Methods

run()

process results using methods from fit process and emit

1.10 SpinorMonitor

Contents:

1.10.1 MainWindow

class SpinorMonitor.**MainWindow**

Main Window for the app, contains the graphs panel and the options panel. Executes main control of all other panels.

Variables

- **expData** – Pandas dataframe where all experiment information is kept
- **run** – Run number for the day
- **path** – Path to data storage folder
- **processThreadPool** – Dictionary of running threads
- **process** – Convenience dictionary to initialize objects
- **ROI** – region of interest
- **running** – Boolean if data collection thread is active
- **index** – keeps track of shot internally
- **image** – ImageWindow widget
- **plots** – DataPlots widget
- **options** – Options widget
- **plot_options** – PlotOptions widget
- **vis_plots** – VisualPlotter widget
- **data_tables** – DataTable widget

- **aux_funcs** – AuxillaryFunctionContainerWidget widget
- **tabs** – QTabWidget, contains other widgets
- **ipy** – QIPythonWidget

Methods

center()
Centers Window

change_state()
start and stop data collection thread

data_process(results_dict)
process the data, including spawn a thread and increment index

data_recieved()
Send message that data was recieved

end()
function to stop listening Thread, writes out expData to csv in smae folder as data printing

finish_thread(ind)
pop the process should destroy it all I think/

get_options()
convenience function to return list of options note that function which recieves params must make deep copy or there will be problems!!

get_roi()
returns region of interest in list :returns: [xstart,xend,ystart,yend,angle] :rtype: list

initUI()
Iniitalize UI and name it. Creat all children widgets and place them in layout

on_fit_name(data)
Triggers the plots.change_key functions with argument data.
Params data name of fit

on_message(data)
Send message to output windows
Parameters data (object) – message to send

set_up_ipy()
setup the ipython console for use with useful functions

start()
Function to start listening thread, connect signals and :var imageThread: IncomingImage object listening for images

to_ipy()
push all variables to Ipython notebook

update_data(results_passed)
function to update plots and push data to ipython notebook

1.11 Optionswidgets

Contents:

1.11.1 FitInfo

class Optionswidgets.**FitInfo** (*params, parent=None*)
custom dialog for fit information

Methods

close ()

parse_params (*tabs*)
populates the tables, row and column determined by run and parameter, so same for all table

1.11.2 Options

class Optionswidgets.**Options** (*parent=None*)
Panel which defines options for fitting and analyzing images

Methods

create_fit_panel ()
create a fit panel

fit_name

get_fit_info ()
popup window which has info of all fits

make_key (*index*)

message

remove_fit_panel ()
remove fit panel

save_params ()
update params

set_current_fit (*fit_name*)

1.11.3 ParameterEntry

class Optionswidgets.**ParameterEntry** (*params, first, parent=None*)
popup box to select parameters

Methods

readout ()
function to return updated Parameters object

1.11.4 PlotOptions

class Optionswidgets.**PlotOptions** (*parent=None*)
Widget for Region of Interest Information and other plot options

Methods

set_roi (*vec*)
Generate roi strings and print coords

1.12 Auxfunctions

Contents:

1.13 Visualplotterwidget

Contents:

1.13.1 ParamEntry

class Visualplotterwidget.**ParamEntry** (*parent=None*)
convenience container widget to hold parameters

Methods

1.13.2 PopPlot

class Visualplotterwidget.**PopPlot** (*mod=None, params=None, do_fit=False, parent=None*)
popup class for plots both static and updating

Variables

- **ax** – matplotlib axis
- **figure** – matplotlib figure
- **canvas** – matplotlib canvas
- **toolbar** – matplotlib navigation toolbar

Parameters

- **mod** (*lmfit.Model*) – lmfit Model object for fitting
- **do_fit** (*Boolean*) – Boolean if fitting should occur
- **params** (*lmfit.Parameters*) – fit parameters

Methods

plot (*x*, *y*, *xl*, *yl*, *title*, *std*)

plot the data with a new fit if `do_fit == True`

Params x x vector of points

Params y y vector of points

Params xl x label

Params yl y label

Params title title of plot

Params std standard deviation of points

update (*x*, *y*, *std=None*)

update the plots call the plot function

Params x x vector of points

Params y y vector of points

Params std standard deviation of points

update_init (*xl*, *yl*, *title*, *ignore*, *start*)

update the parameters to start

Parameters

- **title** (*string*) – title of plot
- **xl** (*string*) – x label
- **yl** (*string*) – y label
- **start** (*int*) – starting index

1.13.3 VisualPlotter

class `Visualplotterwidget.VisualPlotter` (*parent=None*)

Class to choose plotting visually so it is easy. Will also automatically update plots for every shot. Can automatically fit on a single shot or updating shot basis.

Variables

- **message** – pyqtSignal which can be transmitted to main message box
- **plots** – Dictionary to hold all the plots
- **data** – local copy of entire pandas dataframe
- **index** – index of shot
- **start** – start of plot region
- **end** – end of plot region
- **ignore_list** – list of shots to ignore

Fit_models different models to fit too needs to be updated when models added

Methods

add_fitting_widgets ()
function populates stacked box for each type of fit

avg_data ()
average data and transform self.x_data and self.y_data this is a really crappy algorithm, but it does the trick

do_fit ()
do a fit

filter_ignore (data)
filter data, list of indices to remove built list of indices not ignored

ignore_update ()
update the ignore list parse out test

make_title_string ()
make a title string

make_updating_title_string ()
make a title string

message

plot_clicked ()
function called when any plot option is called, sets the start and end values

static_plot ()
create new modal popup static plot

test_fit ()
do a fit on the test plto

test_plot ()
update the test plot

update_plots (df, index)
Update the updating plots whose references are stored in self.plots

Params df pandas dataframe holding data

Params index index of shot

updating_plot ()
create an updating plot and fill it with parameters gathered from current state of widgets

validate (el)
valid to make sure is a single integer or list comprehension and turn list comprehensions into their equivalent definee here since its an object method in QTGUI

var_push (var_list)
add a list of variables to options

verbose_avg (x, y)
average data and transform self.x_data and self.y_data this is a really crappy algorithm, but it does the trick

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

- add_element() (Auxfuncwid-
get.AuxillaryFunctionContainerWidget
method), 7
- add_fitting_widgets() (Visualplotterwidget.VisualPlotter
method), 16
- add_lines() (Dataplots.ImageWindow method), 10
- add_plot() (Ipython.PlotObj method), 8
- AuxillaryFunctionContainerWidget (class in Aux-
funcwidget), 7
- AuxillaryFunctionWidget (class in Auxfuncwidget), 7
- avg_data() (Visualplotterwidget.VisualPlotter method),
16
- BEC_num() (Fitobject.fit_object method), 4
- BEC_num_1() (Fitobject.fit_object method), 4
- bimod2min() (Fitobject.fit_object method), 5
- calculate() (Auxfuncwidget.AuxillaryFunctionWidget
method), 7
- center() (SpinorMonitor.MainWindow method), 12
- change_state() (SpinorMonitor.MainWindow method), 12
- clearTerminal() (Ipython.QIPythonWidget method), 8
- close() (Optionswidgets.FitInfo method), 13
- create_fit_panel() (Optionswidgets.Options method), 13
- create_vecs() (Fitobject.fit_object method), 5
- data_process() (SpinorMonitor.MainWindow method), 12
- data_recieved() (SpinorMonitor.MainWindow method),
12
- DataTable (class in Datatablewidget), 7
- do_fit() (Visualplotterwidget.VisualPlotter method), 16
- end() (SpinorMonitor.MainWindow method), 12
- executeCommand() (Ipython.QIPythonWidget method),
8
- filter_ignore() (Visualplotterwidget.VisualPlotter
method), 16
- FingerTabBarWidget (class in Auxwidgets), 9
- finish_thread() (SpinorMonitor.MainWindow method),
12
- fit_image() (Fitobject.fit_object method), 5
- fit_name (Optionswidgets.Options attribute), 13
- fit_object (class in Fitobject), 3
- FitInfo (class in Optionswidgets), 13
- gauss_2D() (Fitobject.fit_object method), 5
- generate_info_widgets() (Auxfuncwid-
get.AuxillaryFunctionWidget method), 7
- generate_params_widgets() (Auxfuncwid-
get.AuxillaryFunctionWidget method), 7
- get_angled_line() (Fitobject.fit_object method), 6
- get_fit_info() (Optionswidgets.Options method), 13
- get_options() (SpinorMonitor.MainWindow method), 12
- get_params() (Auxfuncwidget.AuxillaryFunctionWidget
method), 7
- get_roi() (SpinorMonitor.MainWindow method), 12
- get_vars() (Ipython.SpinorPlot method), 9
- ignore_update() (Visualplotterwidget.VisualPlotter
method), 16
- ImageWindow (class in Dataplots), 10
- IncomingImage (class in Image), 10
- initUI() (SpinorMonitor.MainWindow method), 12
- line_profile() (Fitobject.fit_object method), 6
- MainWindow (class in SpinorMonitor), 11
- make_key() (Optionswidgets.Options method), 13
- make_title_string() (Visualplotterwidget.VisualPlotter
method), 16
- make_updating_title_string() (Visualplotterwid-
get.VisualPlotter method), 16
- message (Optionswidgets.Options attribute), 13
- message (Visualplotterwidget.VisualPlotter attribute), 16
- multiple_fits() (Fitobject.fit_object method), 6
- newImage() (Image.IncomingImage method), 10
- on_fit_name() (SpinorMonitor.MainWindow method), 12
- on_message() (SpinorMonitor.MainWindow method), 12
- Options (class in Optionswidgets), 13
- output() (Auxwidgets.Textbox method), 9
- paintEvent() (Auxwidgets.FingerTabBarWidget method),
9
- ParamEntry (class in Visualplotterwidget), 14
- ParameterEntry (class in Optionswidgets), 13

parse_params() (Optionswidgets.FitInfo method), 13
partial_TF_2D() (Fitobject.fit_object method), 6
plot() (Visualplotterwidget.PopPlot method), 15
plot_clicked() (Visualplotterwidget.VisualPlotter method), 16
PlotGrid (class in Dataplots), 10
PlotObj (class in Ipython), 8
PlotOptions (class in Optionswidgets), 14
PopPlot (class in Visualplotterwidget), 14
popup() (Dataplots.ImageWindow method), 10
printText() (Ipython.QIPythonWidget method), 8
process_results() (Fitobject.fit_object method), 6
ProcessImage (class in Image), 11
pushVariables() (Ipython.QIPythonWidget method), 8

QIPythonWidget (class in Ipython), 8
QIPythonWidgetContainer (class in Ipython), 8

re_import() (Auxfuncwid-
get.AuxillaryFunctionContainerWidget
method), 7
readout() (Optionswidgets.ParameterEntry method), 13
remove_fit_panel() (Optionswidgets.Options method), 13
run() (Image.IncomingImage method), 10
run() (Image.ProcessImage method), 11

save_params() (Optionswidgets.Options method), 13
set_axis() (Ipython.SpinorPlot method), 9
set_current_fit() (Optionswidgets.Options method), 13
set_roi() (Optionswidgets.PlotOptions method), 14
set_up_ipy() (SpinorMonitor.MainWindow method), 12
setImage() (Dataplots.ImageWindow method), 10
sg2min() (Fitobject.fit_object method), 6
SpinorPlot (class in Ipython), 9
start() (SpinorMonitor.MainWindow method), 12
static_plot() (Visualplotterwidget.VisualPlotter method), 16
stern_gerlach_2D() (Fitobject.fit_object method), 6
subtract_background() (Fitobject.fit_object method), 6

tabSizeHint() (Auxwidgets.FingerTabBarWidget method), 9
test_fit() (Visualplotterwidget.VisualPlotter method), 16
test_plot() (Visualplotterwidget.VisualPlotter method), 16
TextBox (class in Auxwidgets), 9
TF_2D() (Fitobject.fit_object method), 4
Therm_num() (Fitobject.fit_object method), 5
to_ipy() (SpinorMonitor.MainWindow method), 12

update() (Ipython.PlotObj method), 8
update() (Visualplotterwidget.PopPlot method), 15
update_data() (SpinorMonitor.MainWindow method), 12
update_init() (Visualplotterwidget.PopPlot method), 15
update_pandas_table() (Datatablewidget.DataTable method), 7
update_plot() (Ipython.SpinorPlot method), 9
update_plots() (Visualplotterwidget.VisualPlotter method), 16
updatePlot() (Dataplots.ImageWindow method), 10
updating_plot() (Visualplotterwidget.VisualPlotter method), 16

validate() (Visualplotterwidget.VisualPlotter method), 16
var_push() (Visualplotterwidget.VisualPlotter method), 16
verbose_avg() (Visualplotterwidget.VisualPlotter method), 16
VisualPlotter (class in Visualplotterwidget), 15